

Text Document

Generated on: 7/16/2026

CardioTrack CT-200 Blood Pressure Monitor

1. Introduction

The CardioTrack CT-200 is a digital monitor intended for use in measuring blood pressure and pulse rate in adult patient populations.

2. Safety Precautions

WARNING: Do not use this monitor on infants, toddlers, or individuals <-- (ADDED: "toddlers")

who cannot express consent.

WARNING: Avoid operating the device near strong electromagnetic fields, such as microwave ovens or cellular phones.

3. Measurement Guidelines

3.1 Cuff Inflation

Wrap the cuff snugly around your upper left arm, about 1-2 cm above the elbow joint. Press the Start button. The pump will automatically inflate the cuff to 180 mmHg. <-- (MODIFIED: limit raised to 180 mmHg)

3.2 Device Calibration <-- (NEW SECTION ADDED)

The monitor performs automatic self-calibration during each power-on cycle.

Annual calibration verification by an authorized service center is recommended.